

# F-DECALIN PLUS



More Information:  
[www.fluoron.de](http://www.fluoron.de)

## Biocompatible

## Physiologically and chemically inert

## Adjusted density of 1.81 g/cm<sup>3</sup>

F-DECALIN<sup>PLUS</sup> is a mixture of Perfluorodecalin and F6H8<sup>®</sup>. Its reduced density promotes "gentle" reattachment of the retina.

- Intraoperative adjuvant for retinal reattachment.
- The addition of F6H8<sup>®</sup> allows remaining droplets to build a homogenous tamponade medium with the ensuing silicone oil tamponade.
- The reduced density allows tractions to remain visible to the surgeon.

**G-80205 F-DECALIN<sup>PLUS</sup>**  
Vial 5 ml, sterile

**G-80207 F-DECALIN<sup>PLUS</sup>**  
Vial 7 ml, sterile