



FEMTO LASER INSTRUMENTS

**INSTRUMENTS FOR
LASIK AND CATARACT
SURGERIES**



Geuder[®]
Precision made in Germany

INSTRUMENTS FOR CATARACT SURGERY



G-33965

DICK

FEMTOCATARACT EYE SPECULUM

open valves with middle bar,
valve length 18 mm,
total length 80 mm



G-33966

DICK

FEMTOCATARACT SPATULA

double ended,
spatula sharp, 7 mm angled,
2.5 mm recessed,
0.65 mm wide
knife 1 mm angled,
0.9 mm wide,
length 100 mm



G-33967

BREYER

FEMTOCATARACT SPATULA

blunt, 6 mm angled,
tip 1 mm angled, blunt,
length 100 mm




G-33968

HOFFMANN

FEMTOCATARACT NUCLEUS ROTATOR AND CHOPPER

Chopper 1 mm,
bracket 1 mm long,
0.75 mm wide,
length 120 mm






G-32926
AUFFARTH
MICS
DOUBLE CROSS ACTION
CAPSULORHEXIS FORCEPS
 Pure Efficiency Phaco
 flat handle, cross-action,
 for incision size 1.6 mm,
 total length 100 mm



G-15250
SAUTTER
HYDRODISSECTION
CANNULA
 27 gauge / 0.4 mm




G-34245
SAUTTER
SINGLE-USE HYDRO-
DISSECTION CANNULA
 27 gauge / 0.4 x 22 mm
 10 pcs. per box, **sterile**




G-32564
RENTSCH
CAPSULE POLISHING
CURETTE
 angled, 90°

INSTRUMENTS FOR LASIK SURGERY



G-33952
CHANSUE
FEMTO
DOUBLE INSTRUMENT
lenticule spatula, blunt
and hook



G-33950
CHANSUE
FEMTO
DOUBLE INSTRUMENT
lenticule spatula, blunt spoon
and hook



G-33954
GUELL
FEMTO
DOUBLE INSTRUMENT
lenticule spatula,
blunt spoon, rounded tip,
slightly curved
and hook



G-33960
RUPAL SHAH
LENTICULE FORCEPS
curved, serrated tip,
6 mm titanium handle,
blue-colored



G-18914
BURATTO
LASIK CORNEAL FLAP
FORCEPS



G-33961
LENTICULE FORCEPS
bent, 7.5 mm,
with fine grooves



G-33759
FIXATION RING
inside diameter 13 mm,
swivelable, 9 teeth, blunt



G-14010
KUHNT
CORNEAL SCARIFIER
fig. 1



G-33941
KELMAN-MCPHERSON
FRAGMENT FORCEPS
especially designd for
FEMTOLASIK,
angled, 7.5 mm,
with fine grooves




G-33428
PFAEFFL
FLAPORHEXIS FORCEPS
Jaws oval

INSTRUMENTS FOR LASIK SURGERY



G-33429
KNORZ
**FEMTO-LASIK FLAP
FORCEPS**
grasping surface roughened,
curved



G-33940
BLUM
**FEMTOLASIK ROUND
SPOON**
bevel up, semi sharp
angled, with radius




G-33942
BREYER
**FEMTOLASIK
DOUBLE-ENDED
INSTRUMENT**
Blum round spoon,
modified by Breyer
semi-sharp,
hook with one blunt and one
slightly sharp edge




G-33427
PFAEFFL
FEMTOFLAPLIFTER
both tip ends blunt,
curved part blunt



G-33943
PFAEFFL
FEMTOFLAPLIFTER
v-shaped hook with straight
front part



G-33958
PFÄFFL - BREYER
FEMTOFLAPLIFTER
V-spatula with recessed tip



G-33944
KOHNEN
FEMTOLASIK
DOUBLE-ENDED
INSTRUMENT
spatula and blunt hook



G-16189
SPATULA
45°, angled, 0.5 x 0.25 mm
atraumatic



G-33959
BREYER
FEMTOLASIK
HOOK INSTRUMENT
angled

GEUDER AG reserves the right to make changes to technical details in response to recent developments. Geuder does not assume liability for the accuracy of each individual statement.

Illustrations not drawn to scale.



Distributed by:

IQ Medical Pty Ltd

2/86 Mary Street, Unley SA 5061

Phone (08) 8357 8022

Email sales@iqmedical.com.au

Web www.iqmedical.com.au

Geuder[®]
Precision made in Germany

GEUDER AG
Hertzstrasse 4
69126 Heidelberg
Germany

Phone: +49 6221 3066
Fax: +49 6221 303122
info@geuder.de
www.geuder.com