


Only Lenstec offers the precision and predictability of an aspheric IOL made from unique blue-blocking material proven effective for 20 years.



Distributed by:

IQ Medical Pty Ltd

2/86 Mary Street, Unley SA 5061

Phone (08) 8357 8022

Email sales@iqmedical.com.au

Web www.iqmedical.com.au



Just because an IOL looks yellow does not mean it will stay yellow... or that it will absorb light of the correct wavelength to benefit your patient.

...now tested for 20 year stability.*

The Softec HDY is fabricated from a material that has been proven stable, particularly with respect to long term exposure to sunlight, using an accelerated aging study that replicates the equivalent of 20 years of exposure for the average person (based on specific ISO standards for the photostability of IOLs).


So, while the discussion of using blue light absorbers in the fabrication of IOLs has been ongoing since the late 1980s, most researchers would agree that blocking certain wavelengths of visible light is beneficial to the eye... and a material proven to provide stable light transmission properties for 20 years is the obvious choice.

Softec HDY

*High Definition • Micro Incision • Aspheric IOL • Precision Diopters
.25 Diopter Increments from +18.0 to +25.0*

Micro Injection System
via 1.6mm cartridge

.25 Diopter
Increments



Optic Size:	5.75mm
Optic Type:	Equal Conic Bi-aspheric
Length:	12.00mm
Haptic Style:	Modified C
Angulation:	0 Degrees
Construction:	1 Piece
Positioning Holes:	0
Optic Material:	Acrylic (26% Water Content)
A Constant:	118.0
A/C Depth:	5.10mm
Diopter Increments:	
Whole	+5.0 to +36.0
Half	+10.5 to +29.5
0.25	+18.0 to +25.0

U.S. Patent No. 20060279697 A1 December 14, 2006

* see http://www.gclabsite.com/files/publication/pub46_6-9.pdf
The Softec HDY is not currently FDA approved and is not for sale in the US.

Lenstec UK
Lenstec Barbados Inc.
Lenstec House, Lakeside
Calder Island Way, Wakefield
West Yorkshire WF2 7AW
United Kingdom
Tel: +44 (0)1924 382 678
Fax: +44 (0)1924 831 080
Email: lenstecuk@lenstec.com

Lenstec Barbados
Lenstec Barbados Inc.
Airport Commercial Centre
Pilgrim Road
Christ Church BB17092
Barbados
Tel: 246-420-6795
Fax: 246-420-6797
Email: lenstecbarbados@lenstec.com

Lenstec Spain
Lenstec Barbados Inc.
Avda Lluís Companys, 62, office 9
08302 Mataró, Barcelona
Spain
Tel: +34 937570482
Fax: +34 937986329
Email: lenstecsain@lenstec.com

World Headquarters
Lenstec Inc.
1765 Commerce Ave. N.
St. Petersburg, FL 33716
USA
Tel: 727-571-2272
Fax: 727-571-1792
Email: lenstec@lenstec.com

www.lenstec.com



Lenstec is an ISO 13485 Registered company manufacturing CE Marked products.