



# Catalogue Number and Packaging

Provided sterile. 6 packages of 4 MST Capsule Retractors per box.

Description:  
MST Capsule Retractors


Catalogue Number  
MCR 0001



There are 4 retractors in each package, six packages in a box.




The retractors are held in place by a silicone bar and are gripped by slots in the bar.




The Tyvek lid is pulled back and the tray turned upside down. The bar will drop out.



Like this...



The bar is picked up and squeezed as shown to release its grip on the retractors.



Remove the retractors.



Distributed by:  
**IQ Medical Pty Ltd**  
2/86 Mary Street, Unley SA 5061  
Phone (08) 8357 8022 | Email [sales@iqmedical.com.au](mailto:sales@iqmedical.com.au)  
Web [www.iqmedical.com.au](http://www.iqmedical.com.au)