


Maestro and Barrett I/A Catalog Numbers



MST Maestro Handpiece

MHP-0001

Barrett I/A Tips for use with Maestro on Alcon® Systems

Incision	Included Sleeve	Description	Catalog No.
2.6-2.85	 Purple	MST Barrett, Straight Capsule Polish 0.3mm port, w/ 0.7 A1 purple sleeve*	BRT-1028-P
		MST Barrett, Bent Capsule Polish 0.3mm port, w/ 0.7 A1 purple sleeve*	BRT-1029-P
		MST Barrett, Curved Capsule Polish 0.3mm port, w/ 0.7 A1 purple sleeve*	BRT-1030-P
2.2-2.5	 Red	MST Barrett, Straight Capsule Polish 0.3mm port, w/ 0.7 A1 red sleeve*	BRT-1028-R
		MST Barrett, Bent Capsule Polish 0.3mm port, w/ 0.7 A1 red sleeve*	BRT-1029-R
		MST Barrett, Curved Capsule Polish 0.3mm port, w/ 0.7 A1 red sleeve*	BRT-1030-R
1.8-2.4	 Gold	MST Barrett, Straight Capsule Polish 0.3mm port, w/ 0.7 A1 gold sleeve*	BRT-1028-G
		MST Barrett, Bent Capsule Polish 0.3mm port, w/ 0.7 A1 gold sleeve*	BRT-1029-G
		MST Barrett, Curved Capsule Polish 0.3mm port, w/ 0.7 A1 gold sleeve*	BRT-1030-G

Barrett I/A Tips for use with Maestro on AMO® and B + L® Systems

2.6-2.85	 Blue	MST Barrett, Straight Capsule Polish 0.3mm port, w/ 0.7 blue sleeve	BRT-1033-B
		MST Barrett, Bent Capsule Polish 0.3mm port, w/ 0.7 blue sleeve	BRT-1034-B
		MST Barrett, Curved Capsule Polish 0.3mm port, w/ 0.7 blue sleeve	BRT-1035-B
2.4-2.75	 Yellow	MST Barrett, Straight Capsule Polish 0.3mm port, w/ 0.7 yellow sleeve	BRT-1033-Y
		MST Barrett, Bent Capsule Polish 0.3mm port, w/ 0.7 yellow sleeve	BRT-1034-Y
		MST Barrett, Curved Capsule Polish 0.3mm port, w/ 0.7 yellow sleeve	BRT-1035-Y
1.8-2.4	 Gold	MST Barrett, Straight Capsule Polish 0.3mm port, w/ 0.7 gold sleeve*	BRT-1033-G
		MST Barrett, Bent Capsule Polish 0.3mm port, w/ 0.7 gold sleeve*	BRT-1034-G
		MST Barrett, Curved Capsule Polish 0.3mm port, w/ 0.7 gold sleeve*	BRT-1035-G

Maestro™ with Barrett I/A™

Alcon is a trademark of Alcon, Inc.

MST (MicroSurgical Technology) is in no way affiliated with Alcon. All trademarks are property of their respective owners.

The incision size examples contain recommended incision sizes. Surgeons may prefer incision sizes significantly different from those listed depending on incision location, geometry, hardness of the cataract, ocular status, and experience of the surgeon.

*Sleeves work with both MST 0.7 and 0.9 tips



Distributed by | IQ Medical Pty Ltd

2/86 Mary Street, Unley SA 5061

Phone (08) 8357 8022 | Email sales@iqmedical.com.au | Web www.iqmedical.com.au